

Enhancements	V4	V2	V1
Reduced weight (7%)	●		
Capture impedance from standard defibrillation pads while monitoring any lead	●		
Keypad improvements—enhanced tactile feedback and upgrades to increase durability	●		
Enhanced AED analysis algorithm	●		
Restriction of Use of Hazardous Substances (RoHS) compliant	●		
STJ values on 12-lead printout	●	●	
Enhanced data software—capture continuous waveforms from all channels	●	●	
Continuous temperature monitoring	●	●	
AC/DC power	●	●	